

Push[®] med Wrist Brace Splint

Allows wrist to rest without impeding hand function

- Wrist flexibility is limited while the hand is supported in a functional position
- Anatomically-shaped, removable stay along the bottom of the wrist follows the natural contour of the wrist joint and palm
- Brace is easily applied and removed with one hand
- Streamlined contour allows for easy dressing and undressing
- Can be machine washed on gentle cycle at low temperature
- Creates soft, comfortable feel and wicks moisture through the Sympress[™] microfibers



Indications

- Sprains / strains
- Tendinitis of wrist tendon(s)
- Carpal tunnel syndrome
- Overuse symptoms
- Arthritis of the wrist



Measuring / Sizing

Measure around the wrist as shown and choose a size from the chart. Note: If you have a cone shaped forearm you may need to choose a larger size.



Sizing: Right or Left Wrist

SIZE	INCHES	CENTIMETERS
1	5 1/8 - 6	13 - 15
2	6 - 6 5/8	15 - 17
3	6 5/8 - 7 1/2	17 - 19
4	7 1/2 - 8 1/4	19 - 21

Insurance Information

- The Push med Wrist Brace Splint has an assigned L-code of L-3908.
- BraceLab is not a provider. We do not file insurance nor can we provide information on insurance coverage or reimbursement rates. Please contact your insurance for filing information.