



Premium European Arm & Wrist Bracing

Push® med Elbow Brace Epi

Relieves lateral or medial epicondylitis pain

- Delivers localized pressure via skin-compatible viscoelastic polymer pad
- Averts distal slippage, especially during active play
- Avoids blood flow constriction & radial nerve pressure via C-shaped reinforcement



Sizing: One Size

L-code does not apply

“I am an Occupational Therapist and a Certified Hand Therapist.... As someone who was treated literally thousands of these injuries, [the Elbow Brace Epi] is by far the best product to treat medial and lateral epicondylitis conservatively.” — Name

Push® med Elbow Brace

Limits elbow extension, allows full flexion

- Prevents full elbow extension
- Allows one-handed adjustable extension restriction
- Withstands forces at variety of elbow angles due to robust strap configuration
- Allows easy application / removal with upper arm zipper
- Precludes brace slippage via internal silicone strips
- Controls elbow edema with circumferential elastic design



Sizing: Right or Left Elbow

L-code does not apply

SIZE	INCHES	CENTIMETERS
1	8 ⁵ / ₈ - 9 ¹ / ₂	22 - 24
2	9 ¹ / ₂ - 10 ¹ / ₄	14 - 26
3	10 ¹ / ₄ - 11	26 - 28
4	11 - 12 ⁵ / ₈	28 - 32
5	12 ⁵ / ₈ - 14 ¹ / ₈	32 - 36



BraceLab

Embrace Functional Freedom

Available in the US only from www.BraceLab.com
support@bracelab.com · P: 888-235-8221 · F: 888-762-6422



Push® med Wrist Brace

Immobilizes wrist with less bulk, greater comfort

- Provides precise fit via multiple individually shaped volar metal stays (non-removable)
- Sandwiches volar and dorsal stays for secure immobilization via wide elastic band closure
- Avoids pressure to vulnerable bony prominences
- Increases comfort, dressing ease with streamlined, shorter design



Sizing: See chart below



Assigned L-code L3908

Push® med Wrist Brace Splint

Limits flexion, supports wrist in functional position

- Features anatomically shaped, removable and adjustable Thermolyn® stay
- Affords easy application / removal with one hand

Sizing, All Wrist Braces: Right or Left

SIZE	INCHES	CENTIMETERS
1	5 1/8 - 6	13 - 15
2	6 - 6 5/8	15 - 17
3	6 5/8 - 7 1/2	17 - 19
4	7 1/2 - 8 1/4	19 - 21



Assigned L-code L3908

Premium Comfort & Care

Push braces are designed in consultation with medical experts to be anatomically correct and provide optimal support with minimal impact on daily activities. They are made in Europe with premium quality, washable materials for comfort and durability, and have a 90-day guarantee. Most Push products are made of Sympress™, a soft, moisture-wicking microfiber.

As the exclusive US distributor of Push braces, your patients' satisfaction is important to BraceLab. We offer exceptional customer service with rapid response and trained support specialists 9am-5pm EST, Monday-Friday.



Ordering Push® Braces

- Order via website or with phone credit card
- Ask us about volume discounts, or how to establish purchase order/billing account

To help your patients order directly:

- Print them a copy of the [Patient Information Sheet](#) from the relevant product webpage
- Help them determine the correct size
- Let them know they must file their own insurance; BraceLab is unable to assist



BraceLab

Embrace Functional Freedom

Available in the US only from www.BraceLab.com
support@bracelab.com · P: 888-235-8221 · F: 888-762-6422

